



# DTX Studio™ Lab

Version 1.12

## **Computer Guidelines**

# Table of Contents

**Recommended Hardware Specifications** ..... **3**

Operating Systems ..... 3

Devices ..... 3

Requirements ..... 3

# Recommended Hardware Specifications

## Operating Systems

- Windows® 64-bit (Pro and Enterprise edition):
  - Windows® 10
  - Windows® 11

## Devices

- Windows desktop
- Windows notebook

## Requirements

<b>CPU</b>	Quad-core (Intel Core i5 or i7)
<b>RAM</b>	8 GB - 16 GB
<b>Graphics card</b>	NVIDIA 2 GB video RAM
<b>Disk space</b>	25 GB free disk space
<b>Network</b>	Broadband Internet connection with 512 Kbps upload speed.
<b>Monitor</b>	Full HD (1920×1080) or higher



#### Manufacturer

Nobel Biocare AB  
Box 5190, 402 26  
Västra Hamngatan 1, 411 17  
Göteborg, Sweden

Phone: +46 31 81 88 00  
Fax: +46 31 16 31 52

nobelbiocare.com

#### Distributed in Australia by:


Nobel Biocare Australia Pty Ltd  
Level 4/7 Eden Park Drive  
Macquarie Park, NSW 2114 Australia

Phone: +61 1800 804 597

#### Distributed in New Zealand by:

Nobel Biocare New Zealand Ltd  
33 Spartan Road  
Takanini, Auckland, 2105 New Zealand

Phone: +64 0800 441 657



Rx Only

**Canada license exemption** Please note that not all products may have been licensed in accordance with Canadian law.

**For Prescription Use Only. Caution** Federal (United States) law restricts this device to sale by or on the order of a clinician, medical professional or physician.

#### Basic UDI Name

Design software

#### Basic UDI-DI

7332747000002146X